


# Surgical Laser Basics: 5 Tips for the Novice User

**John C. Godbold, Jr, DVM**  
 Stonehaven Veterinary Consulting  
 Jackson, Tennessee



1

▲ A full-thickness, single-pass CO<sub>2</sub> laser incision causes minimal thermal tissue change. *Photo courtesy of Dr. Masahiro Seki.*


The CO<sub>2</sub> laser is the predominant surgical laser used for soft-tissue surgery because of its unique interaction with tissue. CO<sub>2</sub> lasers emit an invisible beam of laser light (10,600 nm) that is highly absorbed by water. The absorbed energy is converted into heat, cells are ruptured, tissue is vaporized, and most of the energy delivered to the tissue is released in a plume of steam and cellular remnants.<sup>1</sup>

The highly focused beam of a CO<sub>2</sub> surgical laser is a versatile and precise instrument that can efficiently vaporize soft tissue while simultaneously sealing small blood vessels, lymphatic vessels, and nerve endings.<sup>2,3</sup>

When properly used, a CO<sub>2</sub> laser creates as little as 100 microns of collateral thermal change in tissue. Within this zone of thermal change lie the benefits of CO<sub>2</sub> laser surgery. Sealed blood vessels reduce hemorrhage and increase visualization of the surgical field, sealed lymphatic vessels reduce postoperative edema and swelling, and sealed nerves lessen postoperative pain.<sup>4-9</sup>

To maximize these benefits and achieve best surgical outcome, the laser power (watts), spot size (diameter of the laser beam contacting the tissue in cm<sup>2</sup>), and delivery mode (continuous, single-pulse, multiple-pulse) must be set appropriately, and the surgeon must use appropriate technique.

**The highly focused beam of a CO<sub>2</sub> surgical laser is a versatile and precise instrument that can efficiently vaporize soft tissue.**



▲ The laser handpiece is held perpendicular to the tissue, with correct distance from the tissue, to ensure a predictable spot size of laser light on the tissue. *Photo courtesy of Dr. Paul J. Sessa.*



▲ Tension is applied to tissue when incising and dissecting. Tissue under tension cuts more quickly, resulting in reduced collateral thermal change and better tissue effect. *Photo courtesy of Dr. William E. Shultz.*

**1 Use a Sharp Laser Knife**  
 Power density (ie, intensity of the laser beam applied to an area of tissue in  $\text{watts/cm}^2$ ) dictates how sharp the knife-like laser beam is. The higher the power density is, the sharper the laser knife and more precise the incision. Power density can be altered by changing either the power or spot size. A higher power density allows for faster hand speed, shorter time of tissue exposure to the laser energy, and reduced tissue damage.<sup>5,10</sup>

**2 Be Quick**  
 The best tissue effect results from exposing tissue to laser energy over the shortest time. Use a combination of power, spot size, delivery mode, and hand speed that allows full-thickness incisions in a single pass (*Figure 1*, previous page).

Do not saw through tissue with multiple partial-thickness cuts. Partial-thickness incisions result in repeated tissue exposure to the laser beam. Each

exposure increases the width of the zone of collateral thermal change in the tissue and contributes to cumulative damage. Sawing through an incision with a CO<sub>2</sub> laser can lead to burned tissue, carbon (char) formation, thermal necrosis, delayed healing, and possible dehiscence.<sup>5,11</sup>

### **3 Be Focused**

Minimize the spot size by directing the laser beam perpendicular to the tissue<sup>8</sup> (**Figure 2**, previous page). Know the correct distance the delivery handpiece should be from the tissue (usually 1-3 mm) to deliver the desired spot size. Avoid backing away from the tissue and defocusing the laser beam unless the goal is to produce a low-temperature, non-evaporative thermal change within the tissue. Defocusing is not appropriate except for contraction of a tissue area for hemostasis of small vessels.<sup>4,10</sup>

### **4 Apply Tension**

Whether incising or dissecting, tissue needs to be under tension while being cut. Tension separates tissue as it is cut, exposing new tissue to the laser beam. Tissue cuts more quickly while under tension, resulting in less collateral thermal tissue change<sup>4,5,8</sup> (**Figure 3**, previous page).

### **5 Use the Apps**

Like apps on a mobile device, CO<sub>2</sub> lasers have software programs. These programs adjust emission of the laser light. The laser light can be emitted continuously, as single pulses, or as a series of pulses. Using these apps can help the surgeon adjust tissue effect:

▶ Continuous delivery: The laser energy

is delivered as long as the foot switch is depressed. Continuous delivery is used for routine incisions and dissection where rapid vaporization of tissue is desired.

- ▶ Single-pulse delivery: One timed pulse of laser energy is delivered when the foot switch is depressed. Single-pulse delivery produces a very brief (10-500 msec) pulse of laser light that can be used to precisely remove very small amounts of tissue. This can be useful for dissection in areas of critical anatomy.
- ▶ Multiple-pulse delivery: A series of pulses of laser light are delivered whenever the foot switch is depressed.<sup>10</sup> Each individual pulse maintains effective power density, but because of the pause between pulses, the rate of tissue removal is slower.<sup>12</sup> Multiple-pulse delivery can be useful to achieve good tissue effect when precision requires a slower hand speed (eg, excision of a meibomian gland tumor).
- ▶ Superpulse: Some CO<sub>2</sub> lasers also give the option of delivering the laser beam in either superpulse (SP) or non-superpulse (NON SP) mode. SP mode produces a series of microseconds-long (<800 μsec) very high-energy (60-80 watts) pulses. The time between each high-energy pulse is longer than the time required for thermal relaxation (cooling) of the tissue. In SP mode, the average power delivered over a period of time is the power selected on the laser control panel.

In devices equipped with SP capability, SP can be turned on or off whether in continuous, single-pulse, or multiple-pulse delivery.

**Tissue cuts more quickly while under tension, resulting in less collateral thermal tissue change.**

Excessive collateral tissue change and char formation occur when tissue absorbs heat faster than it can be released by conduction, convection, circulation, and evaporation. SP can help limit collateral tissue change by allowing increased thermal relaxation of the tissue. Since hemostasis occurs in the zone of collateral tissue change, SP mode can interfere with hemostasis if the zone of reversible collateral tissue change becomes too narrow. Use NON SP mode in more vascular tissue where more hemostasis is needed and SP mode in less vascular tissue.

### Closing Thoughts

The CO<sub>2</sub> surgical laser is an easy tool to adopt in soft-tissue surgery. To do so, beginners need to understand the basic concepts of CO<sub>2</sub> laser-tissue interaction and are encouraged to attend an introductory CO<sub>2</sub> laser workshop. They need

to know their laser's power capacity and how to adjust it, how the laser's control panel affects delivery of the laser beam, and how the laser's delivery handpiece is properly positioned in relationship to the tissue.

With that knowledge, beginners will be able to start using laser parameters recommended in introductory courses and by laser manufacturers to achieve good surgical outcomes. Starting with their first surgery, novices can contribute to those good outcomes by applying the tips mentioned here—be quick, be focused, apply tension, and use the apps! ■■■

Suggested reading and author biography available at [cliniciansbrief.com/Surgical-Laser-Basics](http://cliniciansbrief.com/Surgical-Laser-Basics)

### References

- Berger N, Eeg PH. Fundamentals of laser-tissue interaction. In: *Veterinary Laser Surgery: A Practical Guide*. Ames, IA: Blackwell Publishing; 2006:29-41.
- Fligny I, Wu JS, Samonte BR, Fried MP. Comparative study of laser and scalpel nerve transections. *Lasers Surg Med*. 1992;12(1):43-50.
- Dunié-Mérigot A, Bouvy B, Poncet C. Comparative use of CO<sub>2</sub> laser, diode laser and monopolar electrocautery for resection of the soft palate in dogs with brachycephalic airway obstructive syndrome. *Vet Rec*. 2010;167(18):700-704.
- Holt TL, Mann FA. Soft tissue applications of lasers. *Vet Clin North Am Small Anim Pract*. 2002;32(3):569-575.
- Godbold JC Jr. Laser Surgery, CO<sub>2</sub>. In: Norsworthy GD, ed. *The Feline Patient*, 4th ed. Ames, IA: Blackwell Publishing; 2010:655-659.
- Shelley BA. Use of carbon dioxide laser for perianal and rectal surgery. *Vet Clin North Am Small Anim Pract*. 2002;32(3):621-637.
- Xu Y, Li F, Li Z, et.al. A prospective, randomized controlled trial of circumcision in adult males using the CO<sub>2</sub> laser: modified technique compared with the conventional dorsal-slit technique. *Photomed Laser Surg*. 2013;31(9):422-427.
- Bellows J. Laser use in veterinary dentistry. *Vet Clin North Am Small Anim Pract*. 2002;32(3):673-677.
- Duclos DD, Hargis AM, Hanley PW. Pathogenesis of canine interdigital palmar and plantar comedones and follicular cysts, and their response to laser surgery. *Vet Dermatol*. 2008;19(3):134-41, 2008.
- Peavey GM. Lasers and laser tissue interaction. *Vet Clin North Am Small Anim Pract*. 2002;32(3):517-534.
- Berger N, Eeg PH. Power density and the basic effects of radiant energy on tissue. In: *Veterinary Laser Surgery: A Practical Guide*. Ames, IA: Blackwell Publishing; 2006:19-26.
- Berger N, Eeg PH. Clinical applications of CO<sub>2</sub> laser energy. In: *Veterinary Laser Surgery: A Practical Guide*. Ames, IA: Blackwell Publishing; 2006:139-147.